

# Profosma Manual

---



# Jet Plasma



Brand Name	Profosma
Product Name	Jet Plasma
Model Number	R98
Dimension	21*4*5cm
Main material	PC/ABS/ Ceramics
Maximum Voltage	12kV ~ 13kV
Battery capacity	2400mAh
Charging time	2.5H
Product weight	165g
Application	For beauty use

# **Terms and Conditions**

1. Prior to using the machine, ensure to cleanse your face thoroughly and then dry it. Refrain from applying any cosmetics or other skincare products to prevent potential damage to the machine and to maintain optimal effectiveness.
2. Upon completing the operation, you may proceed to apply suitable skincare products based on any skin concerns that may be present.
3. Once you have finished using the machine, remember to power it off. Utilize only a dry, soft cloth or a non-alcoholic damp cloth to clean the machine.
4. Avoid exposing the machine and its associated accessories to liquids and refrain from placing them in damp environments.
5. Store the machine in a well-ventilated and dry area, keeping it out of reach of children.
6. In case the machine's casing sustains damage or if it inadvertently comes into contact with water, discontinue use immediately and seek guidance from customer service. Do not attempt to disassemble or repair the machine without proper authorization.
7. It is imperative to avoid contact with metal objects; kindly remove all metal items such as rings, earrings, and necklaces before undergoing treatment.
8. Keep the machine away from any moisture, and after skin cleansing, ensure to eliminate any moisture with a dry towel before commencing treatment.

# Treatment Parameters

The below is a reference guide. Can be adjusted according to different operations.

Application	Level	Frequency (Hz)	Mode	Time/week	Treatment
Disinfection	Level1	2/3/4/5/10	Pulse	2-3 times	8 times
Shrink pores	Level 1 /2	5/10	Pulse	3-5 times	6 times
Treat acne	Level 2	4/5/10	Pulse	2 times	5-6 times
Evaporate hair follicle oil	Level 2/ 3	000	Continuous	1 time	4-6 times
Collagen regeneration	Level 2	5/10	Pulse	2 times	6-8 times
Anti-inflammatory acne	Level 1/2	2/3/4/5	Pulse	1 time	3-4 times
Whitening rejuvenation	Level 1/2	000	Continuous	3-5 times	5-7 times

## Instructions

1. Prior to treatment, remove all jewelry and ensure your face and hands are free of moisture post-cleansing.
2. The recommended duration for each treatment session is 5 minutes in automatic mode, with each client's facial treatment concluding within this timeframe. Subsequently, the machine will either pause or shut down for two minutes to prevent overheating due to its high-voltage output. After this brief pause, you can proceed with treating the next client. This process should be repeated after every client is treated.
3. Following the treatment, cleanse your face with clean water and apply gentle skincare products or repair essence.
4. Note: During operation, the high-voltage ionization package will emit plasma gas, which may have a distinct odor. This scent is a normal occurrence. The plasma gas, within its concentrated range, possesses sterilizing properties and does not pose any harm to the body. Please feel at ease when using the machine.
5. Post-treatment, it is typical for the skin to exhibit slight redness, which will naturally diminish over time.

# Precautions

1. Minors or elderly individuals are advised against using this machine and should receive instruction from a trained professional.
2. While in operation, ensure the machine is continuously moved and avoid prolonged exposure in a single skin area to prevent potential injury.
3. Each session should ideally not surpass 5 minutes of use, followed by a mandatory two-minute pause.
4. Refrain from utilizing or storing the machine in a bathroom or moist environment to prevent the risk of electric shock or fire hazards.
5. Avoid contact between the plasma ceramic (Jet disposable) and metal objects.

# Machine Integration



To initiate the operation, briefly press the "Shooting button" to commence working. The screen will showcase "START," preventing the operation of other function options during this time. To pause the work, simply short press the "Shooting button" again. This action will prompt the screen to show "READY," allowing for the operation of other function options at that point.



For Pulse Mode:

- Pulse frequencies available: 2Hz, 3Hz, 4Hz, 5Hz, 10Hz.
- Press the "Mode button" briefly to access the "Pulse Mode."
- Press and hold for 4 seconds to cycle through the 5 different frequencies.
- Press and hold for 4 seconds again to lock in the selected frequency, or let it stay for 4 seconds to automatically secure the choice.

For Continuous Mode:

- Continuous frequency: 000.
- Press the "Mode button" briefly to opt for the "Continuous Mode."

# Spot Plasma



Brand Name	Profosma
Product Name	Spot Plasma
Model Number	R96
Dimension	21*4*5cm
Main material	PC/ABS/ PP
Maximum Voltage	2kV
Battery capacity	2400mAh
Charging time	3.5H
Product weight	165g
Application	For beauty use

# Terms and Conditions

1. Prior to machine usage, apply an anesthetic first.
2. Upon completion of the treatment, apply suitable repair skincare products based on any skin symptoms.
3. Avoid exposing the machine and its associated accessories to liquids or humid environments.
4. Store the machine in a dry and well-ventilated area, keeping it out of reach of children.
5. In the event of machine shell damage or water ingress, discontinue use immediately and seek guidance from customer service. Do not try to disassemble or repair the machine without proper authorization.

## Treatment Parameters

The below is a reference guide. Can be adjusted according to different operations.

### **Pulse & Continuous Mode**

<b>Application</b>	<b>Mode</b>	<b>Level</b>	<b>Frequency (Hz)</b>	<b>needles</b>
Skin spots/ dark acne	Pulse	Level 2/3	10/20/40/80Hz	short/long
Over eyeball	Continuous	Level 2/3	000	long
Small black spot on the arm	Pulse	Level 2	100/200Hz	short/long
Eyelid lifting	Pulse	Level 3	200Hz	short/long
Eyelid lifting	Continuous	Level 3	000	copper/short/long
Fat granule	Continuous	Level 1/2/3	000	short/long
Facial acne	Continuous	Level 1/2	000	short/long
Facial acne	Continuous	Level 1/2	000	short/long
Neck Lift	Pulse	Level 2	2/5/10/20/40Hz	short/copper
Stretch Marks	Pulse	Level 2	2/5/10/20Hz	short/long
Chin	Continuous	Level 1	000	short/copper
Ear Lobe Rejuvenation	Pulse	Level 3	40/80/100Hz	short/long
Hyper-Pigmentation	Pulse	Level 2	40/80Hz	short



## One-shot Mode

Application	Pulse Width	Pulse interval	Level	needles
Fat granule	0.2s/0.4s/0.6s	0.1s/0.2s/0.3s	Level 1/2	short/long
Light Eye bags	0.2s/0.4s/0.6s	1.0s	Level 2/3	short/long
Small black spot on the arm	0.8s/1.0s	0.3s	Level 3	short/long
Upper lip tightening	1.0s	0.1s/0.2s/0.3s	Level 1	short/long/copper
Benign Lesions	1.0s	0.4/0.5/0.6s	Level 2	short
Age Spots/ Sun Damage	0.8s/1.0s	0.1s/0.2s/0.3s	Level 2	short
light scar	0.8s/1.0s	0.1s/0.2s	Level 2/3	short/long
scar	1.0s	0.1s	Level 3	long
Jowl	0.2s	0.4/0.5s	Level 2	long
Wrinkles	0.2s/0.4s	0.1s/0.2s/0.3s	Level 2/3	short/long
Head up removal	0.2s/0.4s/0.6s	0.3s	Level 1/2	short/long

# Machine Integration



To initiate the machine's operation, simply give a quick press to the "Shooting button." The screen will then show "START," disabling the use of other function options. If you need to pause the work, another short press on the "Shooting button" will do the trick. This action will change the screen display to "READY," allowing you to make use of other function options at that moment.

For One-Shot MODE:

- Pulse Width options: 0.2s, 0.4s, 0.6s, 0.8s, 1.0s
- Pulse interval options: OFF, 0.1s, 0.2s, 0.3s, 0.4s, 0.5s, 0.6s, 0.7s, 0.8s, 0.9s, 1.0s
- Selection process: Short press the "Mode button" to opt for the "ONE-SHOT MODE." Continue pressing and holding for 4 seconds to choose from the 5 different pulse widths. Repeat the press and hold for 4 seconds to secure the selection or let it remain for 4 seconds to automatically lock the choice. To select the pulse interval, short press the "Multi-function button."

For Pulse MODE:

- Pulse frequency options: 2Hz, 5Hz, 10Hz, 20Hz, 40Hz, 80Hz, 100Hz, 200Hz
- Selection process: Press the "Mode button" briefly to choose the "PULSE MODE." Maintain the press for 4 seconds to cycle through the 8 different frequencies. Hold again for 4 seconds to lock the selection, or allow it to stay for 4 seconds to automatically secure the choice.

For Continuous MODE:

- Continuous frequency: 000
- To opt for "CONTINUOUS MODE," simply give a short press to the "Mode button."

For LED:

- To control LED function:
  - Begin by pressing and holding the "Multi-function button" for 4 seconds. The icon in the upper left corner will start to flash.
  - Next, provide a short press to turn the LED on or off.
  - For locking the LED status, press and hold the "Multi-function button" for another 4 seconds, or allow it to remain for 4 seconds to automatically lock the LED status.